
# REACT

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**Medtronic**  
Further, Together

# CATETER REACT™

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

### New Devices and Techniques

#### ORIGINAL RESEARCH

### A technical comparison of thrombectomy vacuum aspiration systems

Kurt Yaeger,<sup>1,2</sup> Annabelle Iserson,<sup>1</sup> Paul Singh,<sup>1,3</sup> Jacob Wolf,<sup>2</sup> Ester Vidal,<sup>1</sup> Thomas Oxley,<sup>1,2</sup> Anthony B Costa,<sup>1,2</sup> Johanna T Fifi<sup>1,3</sup>

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# INDICATIONS

**CAUTION:** Federal (USA) law restricts these devices to sale distribution and use by or on order of a physician. Indications, contraindications, warnings and instructions for use for Solitaire™ X Revascularization Device can be viewed at [www.medtronic.com/manuals](http://www.medtronic.com/manuals). Indications, contraindications, warnings and instructions for use for all other products can be found in the product labeling supplied with each device.

## Solitaire™ X Revascularization Device

1. The Solitaire™ X Revascularization Device is indicated for use to restore blood flow in the neurovasculature by removing thrombus for the treatment of acute ischemic stroke to reduce disability in patients with a persistent, proximal anterior circulation, large vessel occlusion, and smaller core infarcts who have first received intravenous tissue plasminogen activator (IV t-PA). Endovascular therapy with the device should be started within 6 hours of symptom onset.
2. The Solitaire™ X Revascularization Device is indicated to restore blood flow by removing thrombus from a large intracranial vessel in patients experiencing ischemic stroke within 8 hours of symptom onset. Patients who are ineligible for IV t-PA or who fail IV t-PA therapy are candidates for treatment.
3. The Solitaire™ X Revascularization Device is indicated for use to restore blood flow in the neurovasculature by removing thrombus for the treatment of acute ischemic stroke to reduce disability in patients with a persistent, proximal anterior circulation, large vessel occlusion of the internal carotid artery (ICA) or middle cerebral artery (MCA)-M1 segments with smaller core infarcts (<70 cc by CTA or MRA, <25 cc by MR-DWI). Endovascular therapy with the device should start within 6–16 hours of time last seen well in patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy.

## Phenom™ Catheter

Phenom™ Catheters are intended for the introduction of interventional devices and infusion of diagnostic or therapeutic agents into the neuro, peripheral, and coronary vasculatures.

## React™ Catheter

The React™ 68 Catheter and React™ 71 Catheter are indicated for the introduction of interventional devices into the peripheral and neuro vasculature.

## Riptide™ Aspiration System

The Riptide™ Aspiration System is intended for use in the revascularization of patients with acute ischemic stroke secondary to intracranial large vessel occlusive disease (within the internal carotid, middle cerebral M1 and M2 segments, basilar, and vertebral arteries) within 8 hours of symptom onset. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy are candidates for treatment.

© 2019 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. UC202001600EN

Medtronic  
9775 Toledo Way  
Irvine, CA 92618  
USA  
Tel: 877.526.7890  
Fax: 763.526.7888